ri-vital®

Gebrauchsanweisung
Operating Manual
Mode d'emploi
Instrucciones de uso

CE 0124

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References:
 Chamberlain, J, Terndrup, T., et al. (1995) Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer. *Annals of Emergency Medicine*, Volume 25, Issue 1, Pages 15-20
 Braim, S., Preston, P., and Smith, R. (1988) Getting a better read on thermometry. *RN*. 1998 Mar; 61(3):57-60
 Munk, P, Woods, S, Leduc, D., et al Canadian Paediatric Society Statement; Temperature Measurement in Paediatrics; *Paediatric Child Health* Vol 5 No 5: July/August, 2000
 Brunner, L. and Suddarth, DI, et al. (1982) *The Lippincott Manual of Nursing Practice*, Third Edition; J. B. Lippincott Company, Philadelphia, PA; 1982; p.1145
 Houdas, Y. and Ring, E. F. J.. *Human Body Temperature, Its Measurement and Regulation*; Plenum Press, NY, 1982; p.81-87.

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IMPORTANT INFORMATION TO CONSIDER BEFORE OPERATING	42
PROPER USE	42
WARNINGS AND PRECAUTIONS	43
OVERVIEW OF FUNCTIONS	47
POWER SUPPLY	50
OPERATING CONTROLS / BASIC SETTINGS	51
PULSOXIMETRY	54
PREDICTIVE THERMOMETER	55
BLOOD PRESSURE MEASUREMENT	57
INSTRUCTIONS FOR CARE	63
TROUBLESHOOTING	64
TECHNICAL DATA	67
SCOPE OF DELIVERY / ACCESSORIES	71
TECHNICAL MEASUREMENT CHECK	72
MANUFACTURER'S DATA REGARDING EMC	74
STANDARDS COMPLIED WITH	77

IMPORTANT INFORMATION TO CONSIDER BEFORE OPERATING

You have purchased a high-quality **Riester** product, which was manufactured in accordance with Directive 93/42 EEC and is continuously subject to strict quality control.

Please read this manual carefully before operating the device and keep it in a safe place.

We will be glad to assist you should you have any questions. You will find our address in this manual. We will provide you the address of our partner upon request.

Please note that all of the devices described in this manual are intended solely for use by appropriately trained personnel.

Please also note that correct and safe functioning of this device is only guaranteed when used with the original **Riester** accessories

PROPER USE

The **ri-vital**[®] is a medical device for measuring arterial blood pressure by means of the oscillometric method and arterial oxygen saturation, as well as optional measurement of oral, axillary and rectal body temperature. The pulse oximeter has been calibrated to display the functional oxygen saturation, which is determined from the length of the pulse wave using two high-precision light-emitting diodes (LED). The device is suitable for measurement in adults, children and infants. It can be used to measure and display the values specified in hospitals, clinics, medical practices and similar facilities. Other use is not intended.

WARNINGS AND PRECAUTIONS

EXPLANATION:

WARNING: Indicates things that could harm the patient or user.

CAUTION: Indicates things that could damage the device.

REMARK: Indicates other important information.

WARNINGS:

- Not suitable for use near flammable anaesthetic gas mixtures of air, oxygen or nitrous oxide.
- Do not use this device in environments where there is danger of explosion.
- Do not use this device in the vicinity of magnetic resonance imaging (MR) scanners. MR scanners can generate currents that lead to burns.
- The device's readings must be linked with clinical symptoms and signs. The device is only intended as a tool for diagnosis.
- If you affix certain sensors with plaster or adhesive tape, do not apply it too tightly in order not to interrupt the patient's blood circulation at the measurement site.
- If you connect this device to another device, check that this combination is functioning correctly and ensure that no danger is posed to the patient or the user by this combination. Anyone who combines devices configures a medical system and is therefore responsible for complying with the provisions standard 60601-1-1.
- Incorrectly placed SpO2 sensors can yield incorrect readings. Please heed the operating manual of the sensors.
- The use of defective sensors can lead to incorrect displays that could harm the patient in the worst case. Therefore, you should not use defective sensors.
- Using sensors with damaged cables can also yield incorrect readings. You should therefore check each cable before use.
- This device can display incorrect readings when near strong electromagnetic fields. These fields are generated by high frequency surgical equipment or mobile phones, for example.
- The device can display incorrect readings when near computer tomography (CT) scanners.
- Only use this device with the included accessories and with **Riester** cuffs and sensors.

- To ensure precise and safe blood pressure measurement, a cuff of suitable size must be used. The cuff's marking must be within the markers after it has been placed
- Do not apply any pressure on the tube or the cuff while taking the measurement. This could result in incorrect measurements or to the measurement being aborted.
- Reading errors can occur during blood pressure measurement in patients who suffer from moderate to severe cardiac arrhythmias.
- If continued blood pressure measurements are to be taken from the same patient, then the cuff must be removed between the measurements in order to check the skin underneath the cuff for possible ischemia, discolouration and/or neuropathy.
- Only operate this device with the GlobTek mains adapter of type GTM21089- 1512-T3, as supplied.
- If continued SpO2 measurements are to be taken from the same patient, then the measurement position must be changed regularly. Please check the measuring position against skin irritations, sufficient perfusion and correct position of the sensor. Please change the position of the sensor at least every 4 hours.
- To avoid measurement errors, the patient must remain still while the measurement is being taken.
- This predictive thermometer must NOT be used without a newly placed probe covers from **Riester**.
- Only use original **Riester** probe covers with this device. The use of other probe covers leads to erroneous temperature readings
- The predictive thermometer and the probe covers are NON-STERILE. Do not use on grazed tissue.
- In order to avoid cross-contamination, only use BLUE thermometer probes to measure ORAL or AXILLARY temperature.
- Only use RED thermometer probes for measuring RECTAL temperature.

- Do not open **ri-vital**® and/or the predictive thermometer. The device does not contain any parts that can be maintained by customers. Opening the device can impair the calibration and will invalidate the warranty.
- Disposal of used probe covers must comply with current medical practices and the local regulations concerning disposal in medical facilities.
- Further information on disposal of the lithium batteries is available in the respective community or from the environmental advisory board responsible.
- Used electrical and electronic devices should not be disposed of with normal household waste, but should be disposed of separately according to national or EU guidelines.
- This device is only intended for use by trained medical persons.

CAUTION:

- The device and the finger sensor should be cleaned with a soft, clean cloth. Please use white spirit or isopropyl alcohol to remove heavy soiling. Other agents, especially acetone-based substances, damage the housing or the display.
- Do not press the key with hard or sharp objects.
- Only accessories intended by **Riester** may be used to operate the device.
- A function tester cannot be used to test the measurement accuracy of the pulse oximeter or its sensors.

REMARKS:

- If the blood pressure cuff is not at heart level, a hydrostatic error of about 0.75 mmHg per centimetre will be produced. If the cuff is lower than heart level, about 0.75 mmHg per centimetre must be subtracted from the displayed value. If the cuff is higher than heart level, about 0.75 mmHg per centimetre must be added.
- If a repair to this device is necessary, please contact **Riester** directly or the dealer from whom you purchased the device.
- We accept back devices that have been delivered by us and will dispose of them in an environmentally friendly manner. Please send the device to our company address (carriage paid by receiver).
- The device sticker can display the following symbols:

CE 0124

CEmark

Do not use in the presence of flammable anaesthetic

Patient connection type BF

Serial number

Used electrical and electronic devices should not be disposed of with normal household waste, but should be disposed of separately according to national or EU guidelines.

Device of protection class II

No alarm function

Date of manufacture

Oral or axillar measuring point

Rectal measuring point

Observe the Operating Manual

Sensor not connected to the patient (displayed in the SpO2 window)

Sensor not connected to the device (displayed in the SpO2 window)

Attach probe cover

Not sterile

Latex free

Single use

DEHP free

OVERVIEW OF FUNCTIONS

4.1 Introduction

ri-vital[®] allows you to measure the patient's oxygen saturation of arterial blood and non-invasive measuring of blood pressure, as well as optional measurement of oral, axillary and rectal body temperature. The TFT colour display shows the functional oxygen saturation (SpO2) of the arterial blood, the pulse, the mean (MAP), systolic (SYS) and diastolic (DIA) levels of the arterial blood pressure, as well as the temperature of the selected region of the body. The sensor used for the SpO2 measurement can be inserted into the front of the device for safekeeping. The optional predictive thermometer is equipped with a removable/replaceable insertion device for the thermometer probe, as well as a storage compartment for the probe cover box. The device has a carrying handle which simultaneously serves as a storage place for the blood pressure cuffs. The device can be optionally equipped with a holder for the **ri-thermo**[®] N ear thermometer.

ri-vital[®]

WARNINGS AND PRECAUTIONS

ri-vital[®]

OVERVIEW OF FUNCTIONS

4.2 Information about body temperature

It is a common misconception that 37 °C is the „normal body temperature“. It is actually the case that 37 °C is the average body temperature. The normal body temperature is in a range that varies with age, gender and measuring point. Furthermore, the body temperature fluctuates over the course of the day. It is usually lower in the morning, higher in the afternoon and goes down a little again in the evening. Other factors that affect body temperature may be: the patient’s particular activity, metabolic rate or medications taken. The normal body temperature also tends to drop with increasing age.

The normal temperatures are listed in the following table according to the age of the patient and measuring point. Temperatures measured at different parts of the body, even if they are measured at the same time, must not be directly compared with one another, as body temperature differs between measuring points.

Temperature measuring points	Normal body temperatures according to patient age			
	0-2 years	3-10 years	11-65 years	>65 years
Ear	97,5° - 100,4°F 36,4° - 38,0°C	97,0° - 100,0°F 36,1° - 37,8°C	96,6° - 99,7°F 35,9° - 37,6°C	96,4° - 99,5°F 35,8° - 37,5°C
Mouth	- -	95,9° - 99,5°F 35,5° - 37,5°C	97,6° - 99,6°F 36,4° - 37,6°C	96,4° - 98,5°F 35,8° - 36,9°C
Heart	97,5° - 100,0°F 36,4° - 37,8°C	97,5° - 100,0°F 36,4° - 37,8°C	98,2° - 100,2°F 36,8° - 37,9°C	96,6° - 98,8°F 35,9° - 37,1°C
Rectum	97,9° - 100,4°F 36,6° - 38,0°C	97,9° - 100,4°F 36,6° - 38,0°C	98,6° - 100,6°F 37,0° - 38,1°C	97,1° - 99,2°F 36,2° - 37,3°C
Armpit	94,5° - 99,1°F 34,7° - 37,3°C	96,6° - 98,0°F 35,9° - 36,7°C	95,3° - 98,4°F 35,2° - 36,9°C	96,0° - 97,4°F 35,6° - 36,3°C

References:
 Chamberlain, J. Terndrup, T., et al. (1995) Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer. *Annals of Emergency Medicine*, Volume 25, Issue 1, Pages 15-20
 Braim, S., Preston, P., and Smith, R. (1988) Getting a better read on thermometry. *RN*. 1998 Mar; 61(3):57-60
 Munk, P. Woods, S. Leduc, D., et al Canadian Paediatric Society Statement; Temperature Measurement in Paediatrics; *Paediatric Child Health* Vol 5 No 5: July/August, 2000
 Brunner, L. and Suddarth, DI, et al. (1982) *The Lippincott Manual of Nursing Practice*, Third Edition; J. B. Lippincott Company, Philadelphia, PA; 1982; p.1145
 Houdas, Y. and Ring, E. F. J.. *Human Body Temperature, Its Measurement and Regulation*; Plenum Press, NY, 1982; p.81-87.

4.3 Pulse signal


A detected pulse beat is indicated with an acoustic signal at the moment of systole; the signal’s pitch changes with the value of the oxygen saturation. This allows you to acoustically track the pulse rate and the course of the oxygen saturation without continuously having to keep an eye of the device’s display. **The pulse signal can be deactivated in the Settings menu (Setup button), see section 6.2.**

4.4 Built-in Li-ion battery

The device’s battery permits mobile, battery-supplied operation. The frequency of the blood pressure measurements affects the operating time without a power connection. With blood pressure measurements taken at 10-minute intervals, the **ri-vital®** will work up to 40 hours in battery-operated mode, **with the optional predictive thermometer approx. 20 hours**. A symbol on the display indicates the battery’s remaining charge. Connect the device to a mains socket with the included mains adaptor to charge it. It takes about 8 hours to recharge a fully empty battery. You can operate the device without problem while it is charging. No memory effect occurs due to the use of Li-ion technology.

4.5 Device front

ri-vital® optionally with Predictive Thermometer or **ri-thermo®** Infrared thermometer

- A Carrying handle
- B Coiler for sensor cable
- C Sensor holder with sensor
- D TFT display with illumination
- E On/Off button
- F Charge indicator
- G Operating buttons
- H Predictive thermometer extension module (optional)
- I Probe cover box (optional)
- J Infrared thermometer extension module (optional)
- K Mains adapter connector
- L Charge indicator

5. POWER SUPPLY

The device can be mains operated or off line with the built-in Li-ion battery. For mains operation, the mains adapter must be connected to a power outlet and be inserted into the left side of the device (K) using the cable. Then the battery is also charged. The yellow charge indicator (F) on the front of the device is illuminated while the device is being charged. If it is not illuminated, the battery is fully charged. For the measurement, the device must additionally be turned on on the front of the device using button (E). If no SpO2 measurement is taken and no blood pressure measurement is started either, then the entire device turns itself off in about 10 minutes. This prevents accidental depletion of the battery.

In battery operation, the device shows the remaining charge with a symbol (L) in the lower part of the display. This display has a large tolerance because various parameters, such as ambient temperature, also affect the battery.

Left side view

Right side view

- M Infrared sensors
- N SpO2 sensor connector
- K Mains adapter connector
- O NIBP cuff tube connector
- P NIBP cuff storage place

6. OPERATING CONTROLS / BASIC SETTINGS

6.1 Basic settings menu

After pressing the On/Off button the **Riester** logo appears on the screen, the Print button must be pressed within 5 seconds to arrive in the Basic settings menu.

Set language

Set time

Set date

Delete memory

Leakage test

Manometer mode

Temperature unit

Time format


Only for technical service personnel!

Only for technical service personnel!

The settings can be highlighted yellow with the ↑ ↓ menu buttons and changed with the Setup button.

To exit the menu, press the ↓ menu button until the End-Back setting is reached and confirm the basic settings with the Setup button.

6.2 Settings menu

After pressing the On/Off button the main screen appears after 5 seconds, you arrive in the Settings menu by pressing the Setup button.

Pulse signal

Activate NIBP mode (adult oder neonate)

Cuff pressure

Exit


The „NEO“ setting must only be used for neonates and infants with an arm circumference of less than 13 cm.

The permanent presetting is only used for easier operation in certain surroundings, e.g. the inflation pressure for patients who are mostly hypertonic can be set higher from the outset.

The settings can be highlighted yellow with the ↑ ↓ menu buttons and changed with the Setup button.

To exit the menu, press the ↓ menu button until the End-Back setting is reached and confirm the settings with the Setup button.

6.3 Main screen

- A) Sensor pulsoximeter Status display Sensor crossed out or sensor with arrow / SpO2 reading Pulse bar display
- B) Pulse and perfusion index display
- C) Systolic blood pressure NIBP error message
- D) Mean arterial blood pressure (MAP) Active NIBP mode (only neonate display)
- E) Diastolic blood pressure
- F) Body temperature display (Only with optional predictive thermometer)
- G) Pulsoximeter plethysmogram
- H) Battery charge status
- I) Charge indicator (is turned off when the battery is fully charged)
- J) Error display (only appears with error message)

7. PULSOXIMETRY

7.1 Sensor connector

The sensor connector (N) is located on the left side of the device. The sensor is recessed to prevent damaging the sensor cable, so that the device housing protects the connector. It is recommended that you leave the sensor on the device and use the coiler (B) and the sensor holder (C) for storage. In order to correctly apply the SpO2 sensor to the patient, please read the corresponding operating manual for this accessory part.

7.2 SpO2 measurement

After the sensor has been affixed to the patient, the segments of the bar indicator (A) will start to pulsate. After a short time the readings (A) (B) will appear on the display. The readings are updated once per second.

A symbol with crossed-through sensor appears in the status display (A) if the sensor is not connected. The sensor symbol with arrow appears if the patient is not connected.

Perfusion index:

Measurement of the perfusion index (PI), this PI allows evaluation of the peripheral perfusion of the tissue measured, from which conclusions can be drawn about the quality of the measurement and the blood circulation.

Level 1	< 0.25 %
Level 2	0.25 % - 0.5 %
Level 3	0.5 % - 1.0 %
Level 4	1.0 % - 2.0 %
Level 5	2.0 % - 4.0 %
Level 6	4.0 % - 8.0 %
Level 7	> 8.0 %

7.3 Pulse signal

With the pulse signal switched on, each detected pulse beat is indicated with a brief acoustic signal. The signal's pitch changes with the value of the oxygen saturation. This pulse signal can also be deactivated in the Settings menu.

Please note: no alarm function

PLEASE NOTE: no alarm function

7.4 Pulse display

If the patient is connected to the pulse oximeter, the pulse reading that is obtained in this manner is always indicated.

If no SpO2 measurement takes place, then after a successful blood pressure measurement the device shows the oscillometrically determined pulse. Optionally, the **ri-vital®** can be ordered with Pearl pulse oximeter technology (a process for pulse amplification and artefact suppression).


8. BODY TEMPERATURE MEASUREMENT WITH THE PREDICTIVE THERMOMETER

8.1 Probe covers – application and removal

- A)** Open the probe cover box by lifting and pulling the tab on the upper corner.
- B)** Insert the box with the probe covers from above into the slot. (In order to avoid infections, never interchange or use opened boxes that have been used with blue or red thermometer probes).
- C)** Remove the thermometer probe from the slide compartment. This automatically switches on the thermometer.
- D)** As a reminder to apply or remove the probe cover, a probe symbol with flashing probe cover is displayed before applying or removing the probe cover.
- E)** Insert the probe tip of the thermometer probe into a probe cover in the box. Press the grip of the thermometer probe until the probe cover noticeably clicks into place.
- F)** Take the relevant temperature reading (oral, axillar or rectal temperature). Please use the respective thermometer probe for the purpose.
The measurement is concluded once the final temperature value is displayed
- G)** Press the upper button on the thermometer probe to eject the used probe cover into a waste container.
- H)** Remove, dispose of and replace the empty box.

8.2 Replacing the thermometer probes / inserts


- A) To protect against infection for measurement of ORAL and AXILLAR temperature, only use BLUE thermometer probes. The RED thermometer probes must only be used for measurement of RECTAL temperature. (In order to avoid infections, never interchange or use opened boxes that have been used with blue or red thermometer probes).
- B) To remove or replace the thermometer probes, hold the insert on both sides with thumb and forefinger of one hand whereby the forefinger is placed on the push-button on the back of the insert.
- C) Press the press-button inwards until the lock is released, then push the insert upwards and pull it off the housing.

- D) Now unlock the BLUE or RED connection piece of the thermometer probe by lifting the tab (1), the connection piece (2) can be removed from the insert.

- E) Before the connection piece of a thermometer probe is reused, the insert should be thoroughly cleaned/disinfected as described and also the thermometer probe.
- F) The connection piece of the thermometer probe, as well as the previously removed connection piece, is attached aligned straight and is pushed into the insert until it noticeably locks in place.


- G) Aligned as shown in the picture on the right with the tip of the slot ahead, the insert is now pushed into the housing until the press-button has noticeably locked into place.
- H) The predictive thermometer is ready for use again.

8.3 Oral and axillary measurement:

- A) For measurements of the ORAL temperature, insert the probe tip into the pocket under the tongue on either side of the tongue frenulum (the vertical tissue fold in the middle of the tongue) and push it towards the back.
- B) Remove the thermometer probe and attach a probe cover. The thermometer switches on automatically.
- C) The symbol indicating the ORAL and AXILLARY mode appears on the display. The symbol for the RECTAL mode cannot be displayed if a BLUE thermometer probe is connected.

8.4 Oral measurement:

- A) For measurements of the ORAL temperature, insert the probe tip into the pocket under the tongue on either side of the tongue frenulum (the vertical tissue fold in the middle of the tongue) and push it towards the back.


NOTE: Accurate readings can only be obtained as shown in one of these two heat pockets. Temperature measurements in other areas of the mouth can lead to inaccurate readings.

- B)** The patient's mouth must be CLOSED.
- C)** Hold the thermometer probe firmly at the respective site until the temperature is displayed. Moving the probe during a measurement can impair the capacity of the thermometer to measure the temperature at the site and increase the duration before a reading is obtained.
Press the upper button on the thermometer probe to eject the used probe cover into a waste container.

8.5 Axillary measurement:

- A)** To measure the axillary temperature, request the patient to raise their arms and then place the tip of the thermometer probe with probe cover directly on the patient's skin.
Press gently to ensure good contact. In order to achieve a highly accurate measurement result, place the tip of the thermometer probe with the probe cover directly on the patient's skin.
- B)** Request the patient to lower their arms and to lie still. Keep the thermometer probe almost parallel to the arm.

- C)** The measurement is concluded once the final temperature value is displayed.
- D)** Press the upper button on the thermometer probe to eject the used probe cover into a waste container.

8.6 Measuring the rectal temperature:

- A)** Make sure that the RED thermometer probe is properly inserted.
- B)** Remove the thermometer probe and attach a probe cover. The thermometer switches on automatically.
- C)** A symbol for the RECTAL mode is displayed if a RED thermometer probe is connected.
- D)** Use a medically safe lubricant, as required.
- E)** Insert the thermometer probe in the patient's rectum. In order to ensure good tissue contact, angle the thermometer probe slightly after insertion. Moving the probe during a measurement can impair the capacity of the thermometer to measure the temperature at the site and increase the duration before a reading is obtained.
- F)** The recommended insertion depth 12-19 mm (0.5-0.75 inch) for adults and 6-13 mm (0.25-0.5 inch) for children.
- G)** If two short acoustic signals may be heard, this means that the device has been switched to the DIRECT mode (the slow mode) for this temperature measurement.
- H)** The measurement is concluded once the final temperature value is displayed.
- I)** Press the upper button on the thermometer probe to eject the used probe cover into a waste container.

8.7 Oral, rectal and axillary measurement:

A) If a „D“ appears on the main screen in the temperature display area, this means that the device has been switched to the direct mode (the slow mode) for this temperature measurement.

The predictive thermometer normally works in the predicting mode in order to deliver fast and accurate temperature measurements. However, if no measuring point is recognised or the temperature does not stabilise, the thermometer automatically switches to the direct mode and then works as a temperature monitor. The predictive thermometer also switches automatically into the direct mode if the ambient temperature is over 35 °C (95 °F).

B) After inserting the thermometer probe in the slide compartment, the measured temperature is shown on the main screen until the thermometer probe is removed again or the **ri-vital®** is switched off.

9. Blood pressure measurement:

9.1 Placing the cuff

Slide the cuff over the left upper arm so that the tube is pointing in the direction of the forearm.

Make sure that the underside of the cuff is situated about 2 to 3 cm higher than the elbow and the rubber tube is sticking out from the cuff on the inside of the arm. The marking must be directly over the artery that runs on the inside of the arm.

9.2 Blood pressure measurement

After turning on the device, press button (E) to start the blood pressure measurement. The device will inflate and determine the systolic, diastolic and mean arterial blood pressure with the next discharge. If necessary, the device will automatically inflate more to be able to measure higher pressures as well.

PLEASE NOTE: Pay attention to the NIBP mode setting.

PLEASE NOTE: no alarm function

Should you need to cancel the measurement, then pressing button (E) once more will result in immediate cancellation of the measurement.

9.4 Measurement reading

After a successful measurement, the three blood pressures are shown on the display. If no pulsoximetric measurement is taken, then the device also shows the pulse determined during the blood pressure measurement.

Any errors are documented with an error number in the field of the systolic reading. Please see section 11 for the meaning of the error numbers. Troubleshooting.

9.4 Instructions for taking correct blood pressure measurements

Arterial blood pressure is a vital parameter that changes exceptionally fast as soon as there is the slightest physical exertion. Therefore, please heed the following instructions for determining a correct and reproducible result.

Avoid eating, smoking as well as any exertion immediately before the measurement. All these factors affect the measurement result. The patient should relax in a calm atmosphere for 5 minutes before the blood pressure measurement.

If the patient is wearing clothing that is tight against his/her upper arm, this item of clothing should be removed.


Every attempt of the patient to support his/her arm can increase the blood pressure. You should therefore ensure the patient has a relaxed posture.

Cuffs that are too narrow or too short yield incorrect readings. Selecting the right cuff is exceptionally important. The cuff size depends on the circumference of the upper arm (measured in the middle). The permissible range is printed on the cuff.

A loose cuff or an air chamber that protrudes from the side yields incorrect readings. There may be no gap between the arm and the cuff because otherwise the measurement result will be affected.

Items of clothing that have been raised upwards on the arm must not constrict the arm. Otherwise, these items of clothing (e.g. a pullover) must be removed. Make sure that the cuff is at heart level. Also ensure that none of the tube connections are kinked. The patient should not speak during the measurement.

9.5 Saving vital data:


By pressing the Mem button all readings displayed on the main screen are saved with date and time.

The data saved can be called up with the Print button, they are displayed in chronological order.

There are 1024 storage spaces for complete data records (i.e. 1024 records each with SpO2, pulse, NIBP data etc.).

When the memory is full, each new storage overwrites the oldest reading.

9.6 Printing the vital data (with optional printer):


The data saved can be called up with the Print button, they are displayed in chronological order. By pressing the Print button again in the Storage view the measurement currently highlighted yellow is printed out.

Another measurement can be marked with the ↑ ↓ menu buttons and printed out with the Print button.

Please switch on the printer!

In order to obtain an optimum printing result, the optional printer must be placed as illustrated in the following pictures.

10. INSTRUCTIONS FOR CARE

10.1 Device:

The device, the finger sensor, the optional predictive thermometer (thermometer probe and the slide compartment) should be cleaned with a soft, clean cloth.

Please use white spirit or isopropyl alcohol to remove heavy soiling. After cleaning the device, shake the thermometer probe to remove excess solution. Thoroughly dry the electric contacts of the thermometer probe.

Other agents, especially acetone-based substances, damage the housing or the display.

WARNING!

The holder of the slide compartment must only be cleaned with a soft, clean alcohol-soaked cloth. No liquids must get inside the device.

PLEASE NOTE!:

This thermometer is supplied non-sterile. Do NOT sterilise this device with ethylene oxide gas or heat, in the autoclave or with any methods that may cause damage to materials.

10.2 Disinfectable bladderless cuffs

Cleaning:

The cuff can be wiped with a damp cloth or, like other cuffs, washed with soap in cold water. Please rinse with clear water in this case. In addition, this cuff can be washed up to 60° in the washing machine. Before using again, it must be ensured that no more liquid is in the cuff or the tube. This can adversely affect the measurement result, as well as damage the blood pressure monitor.

The tubes can be wiped off with a damp cloth or a little ethanol on a cotton cloth.

10.3 Disinfection:

The cuff can be completely immersed in disinfectant solution. Only disinfectants with proven efficacy with regard to national requirements should be used.

Before using again, it must be ensured that no more liquid is in the cuff or the tube. This can adversely affect the measurement result, as well as damage the blood pressure monitor.

PLEASE NOTE!

Disinfectable bladderless cuffs must not be ironed! Never subject the cuffs to intense solar radiation! Do not touch the cuff cover nor the pad with a pointed object as this could cause damage!

11. TROUBLESHOOTING

11.1 Predictive thermometer:

The error message appears in the temperature field.

- A)** If the device does not work perfectly after washing and reconnecting the thermometer probe, rinse all contacts with in deionised water and dry with an air blower.
- B)** Symbols indicate all other possible error and corrective action.
- C)** Device error codes:

E01 – System error during synchronisation

E02 – System error during calibration of the patient thermistor

E03 – System error during calibration of the preheating device thermistor

E04 – System timer error

E05 – Preheating device error

E06 – TFT communication error

E07 – Switch-off error

P01 – Probe configuration error (or no thermometer probe connected)

P02 – Patient thermistor instable for the direct mode or outside the valid range

P03 – Preheating device thermistor instable for the direct mode or outside the valid range.

P04 – Patient thermistor instable for the predictive mode or outside the valid range.

P05 – Preheating device thermistor instable for the predictive mode or outside the valid range.

P06 – The probe tip of the thermometer probe cannot be preheated.

D) The device does not contain any parts that can be maintained by customers.

E) For service work please contact your specialist dealer.

11.2 Blood pressure measurement

The error message appears in the temperature field.

Error messages

Inflation time too long.

Causes: No cuff is connected. The cuff is too loose on the arm. The cuff is not tight.

Solution: Place the cuff again and start a new measurement.

Incorrect cuff.

Causes: Adult cuff in neonate mode. Neonate cuff in adult mode.

Solution: Use the correct cuff.

Arrhythmia

Causes: The device determined implausible readings during the last measurement.

Solution: Take another measurement.

Artefact

Causes: The device has detected too many motion artefacts.

Solution: Take another measurement.

Excessive pressure in the cuff or measurement time too long.

Causes: Pressure was applied on the cuff during the measurement.

Solution: Do not touch the cuff and tube during the measurement.

Leak in the cuff or in the connector cable

Causes: The device recorded an unexpected decrease in pressure during the measurement.

Solution: Use a different cuff and check the cable to the device.

For service work please contact your specialist dealer.

12. TECHNICAL DATA

12.1 Pulsoximeter

Note for clinicians: Information regarding the range of the top wavelengths and maximum optical capacity of the light emitted by the pulsoximeter can be useful to clinicians who are applying photodynamic therapy.

Sensors: All PEARL sensors

Measurement

range: SpO₂ 0 - 100 %

Pulse 30 - 249 1/min

Wavelength for the measurement: 660 nm / 880 nm

Intensity: 3.8 mW maximum

Tolerances: SpO₂ +/- 2 % für 70 - 100 %

(below this unspecified)

Pulse +/- 1% or +/- 1 position

Display: Oxygen saturation

(Updating of the readings once per second), pulse as well as a bar graph for pulsation (display is not proportional to the strength of the pulse wave.) and pulsation index

Calibration: Functional saturation.

12.2 Predictive thermometer

Temperature measurement range: 30 °C to 43 °C (85 °F to 110 °F)

Normal predictive measurement range: 35 °C to 43 °C (95 °F to 109 °F)

Direct mode measurement range: 33 °C to 43 °C (91 °F to 109 °F)

Ambient temperature range: 10 °C to 40 °C (50 °F to 104 °F),

for operation: at 10 - 95% relative humidity, non-condensing.

Typical measurement times (after placement at the measuring point): 8 - 10 seconds (Fever temperatures)

Oral (Standard mode): 6 - 10 seconds

Axillary mode:	8 - 12 seconds
Rectal mode:	10 - 14 seconds
Direct mode (all measuring):	60 - 120 seconds
Typical duration before switching to the direct mode:	
At an ambient temperature over 35 °C (95 °F):	Immediate
If no measuring point is recognised:	60 seconds (after removing the probe from the probe slot)
If the temperature does not stabilise:	70 seconds (after placement at the measuring point)
Water bath accuracy (between 35.5 °C and 42.0 °C):	
All measuring points:	±0,1 °C (±0,2 °F)
Thermometer probe:	Replaceable
Insert:	Removable, replaceable, washable, PC+ABS
Materials:	Probe handle, Flame retardant polycarbonate/polyester mix
Probe shaft:	highly fire-resistant polyester
Probe cable:	Polyurethane sheath with TPE overmould
Probe tip:	Aluminium

ALL MATERIALS OF THE PREDICTIVE THERMOMETER ARE LATEX-FREE.

Body temperature measurement with the optional ear thermometer:

The device can optionally be equipped with a holder for the **ri-thermo®** N thermometer in place of the predictive thermometer.

This holder allows convenient storage of the thermometer and the protective sleeves required for the measuring head.

For more information about using the device and important instructions regarding measuring precision, please refer to the operating manual for the **ri-thermo®** N

12.3 Blood pressure (NIBP)

Disinfectable bladderless cuffs: all cuff sizes from the **Riester** portfolio

Measurement range:
Systole 25 - 280 mmHg
Diastole 10 - 200 mmHg
MAP 15 - 260 mmHg

**Measurement range
for neonates:**
Systole 20 - 155 mmHg
Diastole 5 - 110 mmHg
MAP 10 - 130 mmHg

Measurement precision:
Mean value difference < +/-5mmHg
Standard deviation < 8mmHg

Leakage:
< 6 mmHg /minute.

Overpressure limit:
300 mmHg mode „Adult“
150 mmHg mode „Neo“ (neonates)
Overpressure limit (2nd Level): 330 mmHg mode „Adult“
165 mmHg mode „Neo“ (neonates)

Time per measurement
approx. 25s on average „Adult“ mode max. 170 s
„Neo“ mode max. 80s

12.4 ri-vital® General

Mains operation:	Input: 100VAC-240VAC 0,5A 50/60Hz (Mains adaptor Globtek Typ GTM 21089-1512-T3) Output: 12V VDC 1,25A
Capacity:	11 W max.
Protection class:	II
Moisture class	IPX1
Patient module:	Typ BF for SpO2, predictive thermometer
Battery	Li-Ionen Akku, 7,4 Volt 450 mAh
Operating time:	About 250 measurements
Charging:	About 8 hours
Display:	High-contrast TFT colour display with illumination
Temperatures: Storage:	-5 to 55° C
Operation:	10 to 40 ° C
Humidity:	10...95 % relative humidity, non-condensing
Mechanical stability:	Vibration test pursuant to EN ISO 9919:2005 Drop test pursuant to EN ISO 9919:2005
Dimensions (L x B x H):	203 mm x 167 mm x 232 mm with handle 246 mm x 167 mm x 232 mm with handle + ri-thermo® N 273 mm x 167 mm x 232 mm with handle + predictive thermometer
Weight:	1.500g with ri-thermo® N 1.790g with predictive thermometer
Housing:	Break-proof, easy-to-clean PC + ABS plastic

13 SCOPE OF DELIVERY / ACCESSORIES

13.1 Scope of delivery

ri-vital® device
Mains adapter
Fingerclip sensor
Blood pressure connector tube
Blood pressure cuff, size depending on order number
Operating manual

13.2 Accessories

Finger sensor adults (article number 12661)
patient weight > 20 kg
Finger sensor children (article number 12662)
patient weight > 20 kg
Sleeve sensor neonates (article number 12663)
patient weight > 1 kg
Velcro cuff neonates (article number 156)
Velcro cuff infants (article number 158)
Velcro cuff children (article number 154)
Velcro cuff adults, thin arms (article number 160)
Velcro cuff adults (article number 150)
Velcro cuff adults, large arms (article number 152)
Mobile stand with basket (article number 12660)
Infrared printer extension module (holder and printer) (article number 12667)
Infrared printer (article number 12687)
Rolls of paper for infrared printer package of 20 pcs. (article number 12685)
Extension cable for SpO2 length approx. 1.2 m (article number 12664)
Thermometer probe rectal 2.7 m (article number 12668)
Thermometer probe oral/axillary 2.7 m (article number 12669)

13.4 Probe covers:

ri-thermo® N
Package with 25 disposable probe covers (article number 12700)
Package with 100 disposable probe covers (article number 12701)
Package with 1000 disposable probe covers (article number 12702)
Predictive thermometer
Package with 20 disposable probe covers (article number 12684)
Package with 500 disposable probe covers (article number 12688)

14. TECHNICAL MEASUREMENT CHECK

Germany:

Pursuant to the Ordinance on Operating Medical Devices (Medizinprodukte-Betreiberverordnung, MPBetreibV), a technical check must be conducted in two-year intervals. The check can only be conducted by the manufacturer, the appropriate authorities or persons who comply with the requirements of MPBetreibV §6.

European Union outside of Germany:

The respective legal provisions apply for all EU countries outside of Germany.

Technical measurement check for the predictive thermometer

The probe accuracy can be checked using a calibrated water bath (the probe cover must be attached!).

Thermometers which are out of the error threshold must be checked in Direct Mode. The switch to this mode is done automatically when the probe is placed out of the waterbath for approx. 60s.

Technical measurement check for the blood pressure monitor:

1. Procedure for checking compliance with the error threshold of the reading display of the cuff pressure

1.1 Test equipment

- Solid metal container with a capacity of 500 ml +/- 5%.
- Calibrated control manometer with a maximum permissible error of measurement of 1.0 mmHg [0.1 kPa]
- Pressure generator, e.g. pump ball with discharge valve.
- T pieces and tubes.

1.2 Procedure

- Replace the device's cuff with the metal container (see 1.1a).
Connect the calibrated control manometer (see 1.1b) with the cuff connector (see Fig. 1) by means of a T piece and tubes (see 1.1d).
Connect the additional pump ball (see 1.1c) to the pneumatic system with the help of another T piece.
In order to get to the calibration mode (manometer mode), turn the device off and press the On/Off button, the **Riester** logo now appears on the screen; the Print button must be pressed within 5 seconds in order to arrive in the basic settings.

Set language

Set time

Set date

Delete memory

Leakage test

Manometer mode

Temperature unit

Time format


Only for technical service

Only for technical service

Keep pressing the menu button until Start is highlighted yellow in the manometer mode. After pressing the Setup button the device will be in manometer/calibration mode for 170 s.

- Conduct the test in pressure ranges of at most 50 mmHg between 0 mmHg and the highest value of the display range.
Exceeding of the overpressure can take place at any desired level. The test can be continued after decreasing the overpressure to the desired value.

Press the On/Off button to exit the menu.

- Control manometer (see 1.1b)
- Blood pressure device to be tested
- Metal container (see 1.1a)
- Pump ball with discharge valve (see 1.1c)

Fig. 1: Test equipment to determine the error thresholds of the reading display of the cuff pressure.

15. Manufacturer's data regarding EMC

Electromagnetic emissions:

The device is suitable for use in the electromagnetic environment indicated below. The user must ensure that the device is used in such an environment.

Test	Result / Complied with	Comment
HF-emission	Group 1	The device must radiate electromagnetic energy in order to fulfil its function. The functioning of surrounding devices can be affected by this. IEC 61000-4-2
HF-emission	Class B	The device is suitable for use in all facilities. IEC 61000-4-4 IEC 61000-3-2
harmonic emissions	yes	
Flicker emission	yes	IEC 61000-3-3
Voltage fluctuations	yes	IEC 61000-3-3

Electromagnetic stability (immunity):

The device is suitable for use in the electromagnetic environment indicated below. The user must ensure that the device is used in such an environment.

Test	Test step	Classification	Comment
Electrostatic discharge	+/-8 kV air +/-6 kV contact	+/-8 kV air +/-6 kV contact	The floor should be made of a non-conductive material. The relative humidity should be at least 30%. IEC 61000-4-2
Voltage peaks	+/-2 kV mains +/-1 kV input	+/-2 kV mains +/-1 kV input	The mains power should correspond to the power in a typical business environment.

	Output	Output	IEC 61000-4-4
Current surge	+/- 1kV difference +/- 2kV common mode IEC 61000-4-5	+/-1kV difference +/- 2kV common mode	The mains power should correspond to the power in a typical hospital or business environment.
Voltage drops	←5% Un(→ 95% drop in Ut)	see test step	The mains power should correspond to the power in a typical hospital environment. If uninterrupted operation must absolutely be guaranteed, then an uninterrupted power supply must be used.
Test	Test step	Classification	Comment
Voltage drops	40% Un (60% drop in Ut) with five full changes.	see test step	The mains power should correspond to the power in a typical hospital environment. If uninterrupted operation must absolutely be guaranteed, then an uninterrupted power supply must be used.
Voltage drops	70% Un (30% drop in Ut) with 25 full changes	see test step	The mains power should correspond to the power in a typical hospital environment. If uninterrupted operation must absolutely be guaranteed, then an uninterrupted power supply must be used.
Voltage drops	←5% Un (60% drop in Ut) in five seconds	see test step	The mains power should correspond to the power in a typical hospital environment. If uninterrupted operation must absolutely be guaranteed, then an uninterrupted power supply must be used. IEC 61000-4-11
Power frequency	3A/m	3A/m	In the event of operating problems, the device must be placed at a greater distance to the source of the magnetic fields. IEC 61000-4-8

NOTE: Un is the mains voltage for the test device

Electromagnetic stability (immunity), portable RF devices

The distance between mobile communication devices and all parts of the device including all cables should not fall below the distance calculated with the corresponding equation for the frequency of the communication device.

Test	IEC 60601-1-2	Class	Distance
Emitted energy	3V/m	3V/m	Distance = $1.2 \cdot \sqrt{\text{power}}$ for 80.. 800 MHz
	3V/m	3V/m	Abstand = $2.3 \cdot \sqrt{\text{power}}$ for 800 MHz to 2.5 GHz IEC 61000-4-3
Grid-bound energy	3 VRMS	3 VRMS	square root (power) for 150kHz to 80MHz IEC 61000-4-6

The field strengths of stationary senders, such as base stations for cellular phones, of land mobile radio, amateur radio, AM and FM radio transmissions and TV broadcasts cannot be predicted.

A site survey should be conducted for assessing the electromagnetic environment of stationary RF transmitter. If the sender's measured field strength at the device's usage site exceeds the above-mentioned values, then the device's function should be monitored.

If monitoring shows that the device is not functioning properly, corresponding measures must be taken. For example, the device must be realigned or placed at a different location.

Recommended distances

The following minimum distances are recommended between the device and portable HF communication equipment (e.g. mobile phones):

HF frequency	26Mhz - 80MHz	80MHz - 800 MHz	800MHz - 2.5 GHz
Equation	$d=1.2 \cdot \sqrt{\text{power (P)}}$	$d=1.2 \cdot \sqrt{\text{power (P)}}$	$d=2.3 \cdot \sqrt{\text{power (P)}}$
Nominal power (W)	Distance (m)	Distance (m)	Distance (m)
0.1	0.40	0.40	0.80
1	1.20	1.20	1.20
10	3.80	3.80	7.30
50	8.50	8.50	16.30

P is the nominal power of the transmitter in watts, d is the distance in metres

Example:

At a transmission frequency of 700 MHz and power of 10 watts, a distance of $d = 1.2 \cdot \sqrt{10} = 3.79$ m is recommended.

Standards complied with:

Harmonic emissions	IEC 61000-3-2
Flicker emissions	IEC 61000-3-3
Immunity from table electrostatic devices discharge	IEC 61000-4-2, level 3,
Immunity form electromagnetic fields caused by high-frequency radiation	IEC 61000-4-3 at 3V/m
Immunity from sudden voltages peaks/surges	IEC 61000-4-5, level 3
Current surge immunity	IEC 61000-4-5, level 3
Electromagnetic compatibility caused by conduction	IEC 61000-4-6 at 3V/m
Power frequency magnetic fields	IEC 61000-4-8 at 3 V/m
Operation durin mains voltage flucations	IEC 61000-4-11, table 7

Warning:

Using accessories other than those specified can lead to increased emissions and/or decreased stability of the device.

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Références :
 Chamberlain, J, Terndrup, T., et al. (1995) Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer. *Annals of Emergency Medicine*, Volume 25, Issue 1, Pages 15-20
 Braim, S., Preston, P., and Smith, R. (1988) Getting a better read on thermometry. *RN*. 1998 Mar; 61(3):57-60
 Munk, P, Woods, S, Leduc, D., et al Canadian Paediatric Society Statement; Temperature Measurement in Paediatrics; *Paediatric Child Health* Vol 5 No 5: July/August, 2000
 Brunner, L. and Suddarth, DI, et al. (1982) *The Lippincott Manual of Nursing Practice*, Third Edition; J. B. Lippincott Company, Philadelphia, PA; 1982; p.1145
 Houdas, Y. and Ring, E. F. J.. *Human Body Temperature, Its Measurement and Regulation*; Plenum Press, NY, 1982; p.81-87.

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Referencias:

Chamberlain, J, Terndrup, T., et al. (1995) Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer. *Annals of Emergency Medicine*, Volume 25, Issue 1, Pages 15-20
 Braim, S., Preston, P., and Smith, R. (1988) Getting a better read on thermometry. *RN*. 1998 Mar; 61(3):57-60
 Munk, P, Woods, S, Leduc, D., et al Canadian Paediatric Society Statement; Temperature Measurement in Paediatrics; *Paediatric Child Health* Vol 5 No 5: July/August, 2000
 Brunner, L. and Suddarth, DI, et al. (1982) *The Lippincott Manual of Nursing Practice*, Third Edition; J. B. Lippincott Company, Philadelphia, PA; 1982; p.1145
 Houdas, Y. and Ring, E. F. J.. *Human Body Temperature, Its Measurement and Regulation*; Plenum Press, NY, 1982; p.81-87.

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Rudolf Riester GmbH | P.O. Box 35 | Bruckstraße 31 | DE - 72417 Jungingen | Germany